

# Total Cholesterol and Cholesterol Species Determination

Wei Zou (✉ [wzou@dtsc.ca.gov](mailto:wzou@dtsc.ca.gov))

Zou's lab, DTSC, California State Government

---

## Method Article

**Keywords:** cholesterol

**Posted Date:** July 28th, 2011

**DOI:** <https://doi.org/10.1038/protex.2011.245>

**License:**  This work is licensed under a Creative Commons Attribution 4.0 International License.

[Read Full License](#)

---

# Abstract

Total cholesterol and cholesterol species analysis are critical in cardiovascular disease research. The protocol shows procedures that can be used for analysing tissue lipid extracts, lymph, bile or serum.

## Reagents

- Free cholesterol standard solution \((1 \text{ mg/mL in ethanol})
- Cholesterol palmitate standard solution \((1 \text{ mg/mL in chloroform})\): Add 106.383 mg of cholesterol palmitate \((94\%, \text{Sigma})\) into a volumetric flask, top with chloroform.
- 33% KOH
- Ethanol
- Hexane
- O-phthalaldehyde in glacial acetic acid \((50 \text{ mg/100mL})\) Prepared freshly on the day of use
- Concentrated  $\text{H}_2\text{SO}_4$
- Acetone-absolute ethanol \((1:1, \text{v/v})
- Acetone-anhydrous ether \((1:2, \text{v/v})
- Digitonin solution \((0.5\%)\): 100 mg digitonin/20 mL of 50 % ethanol \((10 \text{ mL of D.I. water} + 10 \text{ mL ethanol})\) Make the digitonin solution just before analysis\!
- Glacial acetic acid \((10\% \text{ in water})\): 1 mL of acetic acid in 9 mL of D.I. water

## Procedure

- I. Total cholesterol analysis
1. For each new session:
    - Prepare fresh ethanol, hexane, D.I. water, and o-phthalaldehyde solution.
    - Prepare each set of 16x100 and 16x125 mm test tubes. Check the mouth of tubes.
    - Check the accuracy of automatic pipettes and syringes.
  2. Set up tubes for a standard curve \((duplicate)\):
    - Add separately 10  $\mu\text{L}$ , 20  $\mu\text{L}$ , 30  $\mu\text{L}$ , 40  $\mu\text{L}$  and 50  $\mu\text{L}$  of cholesterol palmitate standard solution to each of the tubes.
    - Evaporate the solvent under the  $\text{N}_2$  in a 40 C water bath.
    - 3. Add samples into a 16x100 mm PTFE capped tube:
    - For tissue lipid extract: use a syringe to add 100  $\mu\text{L}$  of lipid extract and evaporate the solvent under the  $\text{N}_2$  in a 40 C water bath.
    - For lymph/bile/serum: vortex to homogenize the sample, use a 100  $\mu\text{L}$  automatic pipette to add 100  $\mu\text{L}$  of sample, and wash the tip with the same amount of D.I. water.
    - 4. Use an automatic pipette to add 300  $\mu\text{L}$  of KOH into the tube.
    - 5. Use a 5 mL dispenser to add 3 mL of ethanol into the tube, cap it, and vortex it vigorously.
    - 6. Saponify for 15 min in a 60 C water bath and let it cool.
    - 7. Prepare o-phthalaldehyde solution during saponification.
    - 8. Use a 5 mL dispenser to add 5 mL of hexane and vortex vigorously.
    - 9. Use a 5 mL dispenser to add 1.5 mL of D.I. water and vortex vigorously.
    - 10. Transfer the upper layer \((\text{lipid extracts})\) into a 16x125 mm disposable tube w/o cap, wash with 1.5 mL of hexane, combine the extracts and washings, and evaporate the solvent under the  $\text{N}_2$  in a 40 C water bath.
    - 11. Use a 5 mL dispenser to add 3 mL of the o-phthalaldehyde solution into the tube and mix well.
    - 12. Add 1.5 mL of concentrated  $\text{H}_2\text{SO}_4$  by running down the inside wall of the tube and vortex well immediately.
    - 13. Determine the absorbance at 550 nm between 10 and 90 min.
- Wash the cuvettes with D.I. water and methanol to make sure that they are clean and dry for the next sample.
- II. Free cholesterol analysis
1. For each new session:
    - Prepare fresh digitonin solution.
    - Prepare 15 mL centrifuge tubes.
    - Check the accuracy of automatic pipettes and syringes.
  2. Add 1 mL of acetone-absolute ethanol \((1:1)\) into a 15 mL plastic centrifuge tube. Cap the tube.
  3. Add samples into the tube:
  - For tissue lipid extract: use a syringe to add 100  $\mu\text{L}$  of lipid extract.
  - For lymph/bile/serum: vortex to homogenize the sample, use a 100  $\mu\text{L}$  automatic pipette to add 100  $\mu\text{L}$  of sample, wipe the outside of the

pipette tip with Kim-wipers, and wash the tip with acetone-absolute ethanol \(1:1). 4. Use an automatic pipette to add 20  $\mu$ L of glacial acetic acid solution \((10\%) into the tube. 5. Use an automatic pipette to add 500  $\mu$ L of freshly made digitonin solution \((10\%) into the tube. Cap the tube and vortex vigorously. 6. Let the tube stand at room temperature for 1 h. 7. Centrifuge the tube at 1000  $\times$  g for 20 min. 8. Use a Pasteur pipette to remove the supernatant \((cholesterol esters)\) carefully without disturbing the precipitate. 9. Add 1 mL of acetone-ether \((1:2) to the precipitate \((free cholesterol)\). Cap the tube and vortex. 10. Centrifuge the tube at 1000  $\times$  g for 20 min. 11. Use a Pasteur pipette to remove the supernatant \((cholesterol esters)\) carefully without disturbing the precipitate. 12. Dissolve the precipitate in 1 mL of ethanol. Vortex. 13. Transfer the mixture into a scintillation counting vial. Wash the centrifuge tube with 1 mL of ethanol. 14. Evaporate the solvent under  $N_2$  in a 40 C water bath. 15. Redissolve the free cholesterol in 100  $\mu$ L of ethanol. 16. Add 10 mL of scintillation cocktail and count  $^{14}C$  in a liquid scintillation counter.

III. Principles The system, LS6500, is designed for counting  $^{14}C$ ,  $^3H$ ,  $^{32}P$ , and other soft-beta emitters. Radioactive samples measured with an H-number are subjected to quench correction to determine actual sample activity, so called DPM, disintegrations per minute. The standard curve is plotted with the data from quenched standards with known activity containing the same type of isotope as in the sample and H-number are measured for each standard.

IV. Calculations • % E \((efficiency) = \frac{(\text{CPM-background})/\text{DPM}}{\text{DPM}} \times 100\) Where, % E is calculated from correction curve, CPM is obtained in sample, DPM is  $3.92 \times 10^4$  \([= (2.2 \times 10^6 \text{ DPM}/\mu\text{Ci}) \times (0.018 \mu\text{Ci})]\)

• Based on quench correction curve, %E is obtained to calculate DPM for any radioactive samples containing  $^{14}C$ .

• External standard curve: In order to plot standard curve for  $^{14}C$ , quenched standards with known activity were measured. [See figure in Figures section.](#)

## References

- Zou, W., Noh, S.K., Owen, K.O., & Koo, S.I. \((2005)\) Dietary carnitine enhances the lymphatic absorption of fat and a-tocopherol in ovariectomized rats. *J. Nutr.* 135: 753-756.  
["http://jn.nutrition.org/content/135/4/753.full":http://jn.nutrition.org/content/135/4/753.full](http://jn.nutrition.org/content/135/4/753.full) Rudel, L.L. & Morris, M.D. \((1973)\) Determination of cholesterol using o-phthalaldehyde. *J. Lipid Res.* 14: 364-366.
- Sperry, W.M. and Webb, M. \((1950)\) A revision of the Scholtenheimer-Sperry method for cholesterol determination. *J. Biolo. Chem.* 187: 97-100. Wang, C.H., Willis, D.L. & Loveland, W.D. \((1975)\) Radiotracer methodology in the biological environmental and physical sciences. Prentice-Hall, Inc. p. 181-199.

## Acknowledgements

Thanks to Koo, S.I. for mentoring and financial supporting efforts.

## Figures

CPM	H#	%E	DPM
37111	22	94.569	39200
36625	62	93.329	39200
35137	107	89.533	39200
31848	160	81.365	39200
30850	199	78.597	39200

**Figure 1**

Table 1 Results table