

Purification of His-tagged membrane proteins from detergent-solubilized membranes

Dimitrios Fotiadis (✉ dimitrios.fotiadis@ibmm.unibe.ch)

Institute of Biochemistry and Molecular Medicine, and Swiss National Centre of Competence in Research (NCCR) TransCure, University of Bern, CH-3012 Bern, Switzerland

Daniel Harder

Institute of Biochemistry and Molecular Medicine, and Swiss National Centre of Competence in Research (NCCR) TransCure, University of Bern, CH-3012 Bern, Switzerland

Method Article

Keywords: protein purification, membrane protein

Posted Date: August 7th, 2012

DOI: <https://doi.org/10.1038/protex.2012.034>

License: © ⓘ This work is licensed under a Creative Commons Attribution 4.0 International License.

[Read Full License](#)

Abstract

This protocol describes a method to purify a His-tagged membrane protein from detergent-solubilized *Escherichia coli* (*E. coli*) membranes. Solubilized membranes are incubated with Ni-NTA that binds the His-tagged protein. After washing, His-tagged proteins are eluted with histidine. When starting with solubilized membranes or membrane pellets before detergent solubilization, protein purification takes half or one day, respectively.

Reagents

Tris, NaCl, n-dodecyl- β -D-maltoside (DDM) or another detergent of choice, histidine, glycerol, NaN₃, Ni-NTA Superflow beads (Qiagen), Wizard Midicolumn (Promega),

Equipment

SDS-PAGE facility

Procedure

1. Prepare solubilized membranes according to our protocol: Preparation of detergent-solubilized membranes from *Escherichia coli* or start with frozen membranes from 1 l culture by solubilizing in 1% DDM, 20 mM Tris-HCl, pH 8.0, 300 mM NaCl, 10% glycerol, 0.01% NaN₃ for 2 h at 4°C under gentle agitation (final volume: 7 ml) and subsequent ultracentrifugation at 100,000g (4°C, 50 min).
2. Dilute the supernatant 2-fold with Washing Buffer (20 mM Tris-HCl, pH 8, 300 mM NaCl, 0.04% DDM, 5 mM histidine, 10% glycerol, 0.01% NaN₃) and bind for 2 h at 4 °C to Ni-NTA Superflow beads (Qiagen; 0.5 ml bed volume).
3. Load the beads onto a gravity flow column (Wizard Midicolumn, Promega), wash with washing buffer (20 ml), and elute with the same buffer containing 200 mM histidine.
4. Collect fractions of the eluate (e.g. 250 μ l fractions) and identify the fractions containing the purified target protein by SDS-PAGE.
5. Purified protein can be used e.g. for SPA-binding experiments after removing the histidine by desalting columns.

Timing

When starting with solubilized membranes or membrane pellets before detergent solubilization, protein purification takes half or one day, respectively.

References

Casagrande, F. et al. Projection structure of a member of the amino acid/polyamine/organocation transporter superfamily. *J. Biol. Chem.* 283, 33240-33248 (2008).