

# The VetMAX™ *M. tuberculosis* Complex PCR Kit detects MTBC DNA in antemortem and postmortem samples from white rhinoceros (*Ceratotherium simum*), African elephants (*Loxodonta africana*) and African buffaloes (*Syncerus caffer*)

Wynand Johan Goosen (✉ [wjgoosen@sun.ac.za](mailto:wjgoosen@sun.ac.za))

Stellenbosch University Faculty of Medicine and Health Sciences <https://orcid.org/0000-0001-6614-9084>

Tanya Jane Kerr

Stellenbosch University Faculty of Medicine and Health Sciences

Léanie Kleynhans

Stellenbosch University Faculty of Medicine and Health Sciences

Peter Buss

South African National Parks

David Cooper

Enzemvelo KZN Wildlife

Robin Mark Warren

Stellenbosch University Faculty of Medicine and Health Sciences

Paul D van Helden

Stellenbosch University Faculty of Medicine and Health Sciences

Björn Schröder

Thermo Fisher Scientific

Sven David Charles Parsons

Stellenbosch University Faculty of Medicine and Health Sciences

Michele Ann Miller

Stellenbosch University Faculty of Medicine and Health Sciences

---

## Methodology article

**Keywords:** African buffaloes, African elephants, Mycobacterium bovis, Mycobacterium tuberculosis, PCR, VetMAX™ MTBC Detection Kit qPCR, white rhinoceros

**Posted Date:** February 18th, 2020

**DOI:** <https://doi.org/10.21203/rs.2.23897/v1>

**License:** © ⓘ This work is licensed under a Creative Commons Attribution 4.0 International License. [Read Full License](#)

---

**Version of Record:** A version of this preprint was published at BMC Veterinary Research on June 29th, 2020. See the published version at <https://doi.org/10.1186/s12917-020-02438-9>.

## Abstract

**Background:** Bovine tuberculosis and tuberculosis are chronic infectious diseases caused by the *Mycobacterium tuberculosis* complex members, *Mycobacterium bovis* and *Mycobacterium tuberculosis*, respectively. Infection with *M. bovis* and *M. tuberculosis* have significant implications for wildlife species management, public health, veterinary disease control, and conservation endeavours.

**Results:** Here we describe the first use of the VetMAX™ *Mycobacterium tuberculosis* complex (MTBC) DNA quantitative real-time polymerase chain reaction (qPCR) detection kit for wildlife samples. DNA was extracted from tissues harvested from 48 African buffaloes and MTBC DNA was detected (test-positive) in all 26 *M. bovis* culture-confirmed animals with an additional 12 PCR-positive results in culture-negative buffaloes (originating from an exposed population). Of six MTBC-infected African rhinoceros tested, MTBC DNA was detected in antemortem and postmortem samples from five animals. The PCR was also able to detect MTBC DNA in samples from two African elephants confirmed to have *M. bovis* and *M. tuberculosis* infections (one each). Culture-confirmed uninfected rhinoceros and elephants' samples tested negative in the PCR assay.

**Conclusions:** These results suggest this new detection kit is a sensitive screening test for the detection of MTBC-infected African buffaloes, African elephants and white rhinoceros.

## Background

Bovine tuberculosis (bTB) and tuberculosis (TB) are chronic infectious diseases caused by the *Mycobacterium tuberculosis* complex (MTBC) members, *Mycobacterium bovis* (*M. bovis*) and *Mycobacterium tuberculosis* (*M. tuberculosis*), respectively. Many high TB burden countries are dependent on animal-related industries such as agriculture and tourism to support their economy. With encroachment of human settlements into land previously used for agriculture or natural habitats, there are increased opportunities for disease transmission at the interface between animals and people (1). In South Africa, African buffaloes (*Syncerus caffer*) are believed to be a key *M. bovis* wildlife maintenance host, with reported spill-over to other wildlife species and cattle (2,3). In addition, cases of bTB have also been reported in people in African countries, especially those in contact with animals (4). Cases of *M. tuberculosis* infection also occur in domestic cattle in rural areas of Africa as well as captive wildlife and pets (5–8).

In 2016, bTB was discovered in several free-ranging white rhinoceroses (*Ceratotherium simum*) and a single black rhinoceros (*Diceros bicornis*) as well as disease associated with *M. tuberculosis* in an African elephant bull in the Kruger National Park (KNP) (7–9). In South Africa, the largest populations of rhinoceros and elephants reside in KNP, a bTB endemic area. Therefore, infection with *M. bovis* and *M. tuberculosis* have significant implications for species management, public health, veterinary disease control, and conservation endeavours.

Animals infected with MTBC can remain asymptomatic for months or years. Therefore, strategies for management and control require early detection of infection, typically through the measurement of antigen-specific cell-mediated immune (CMI) responses, rather than basing decisions on the presence of disease. Unfortunately, species-specific *in vitro* diagnostic assays for detection of CMI responses to *M. bovis* and *M. tuberculosis* are either very limited, as in the case of white and black rhinoceros (10), or non-existent yet for African elephants. Therefore, antemortem TB diagnosis in these species still relies on mycobacterial culture of respiratory tract samples such as bronchoalveolar lavage fluid (BALF) from rhinoceros or BALF and trunk washes (TW) from African elephants. Because mycobacterial culture is regarded as an imperfect gold standard (11,12), detection of mycobacterial DNA in animal samples has been used as an alternative or ancillary diagnostic method (13). In humans, quantitative real-time polymerase chain reaction (qPCR) assays have become a common molecular tool for rapid detection of TB patients, and now can provide a definitive diagnosis (14,15).

In this study, we describe the use of the VetMAX™ MTBC qPCR Kit for the detection of the insertion element IS6110 from stored tissue and respiratory tract samples previously collected from African buffaloes, white rhinoceros and African elephants with known MTBC-infection status.

## Results

When tissue samples from the randomly selected buffaloes (n = 48) from endemic areas were tested with the VetMAX™ MTBC PCR kit, 38 (79.2%) tested positive (MTBC DNA detected) and 10 (20.8%) tested negative (MTBC DNA not detected). All 26 culture-confirmed *M. bovis*-infected buffaloes were PCR positive (Table 1). Of the 22 culture-negative buffaloes that had been culled based on a positive skin and/or cytokine assay result, 12 were identified as PCR positive (Table 1 and Supplementary Table 1).

Unlike the buffalo, various samples were opportunistically collected from rhinoceros that did not have any previous testing for TB, although they were in a bTB endemic area (KNP). The BALF samples collected from four rhinoceros (animal # 4, 5, 6 and 7) suspected of being *M. bovis*-infected, but unconfirmed by culture, were all PCR negative. Of the six rhinoceros with postmortem samples tested by PCR, MTBC DNA were detected in five animals (Table 1). The 11 tissue samples (lymph nodes and lung) that were *M. bovis* culture-positive from these five rhinoceros were all positive by PCR. An additional two lymph node samples from the same group of *M. bovis*-infected rhinoceros were PCR positive but culture negative. An additional ten tissue samples, collected from four rhinoceros previously identified as *M. bovis*-infected by VetMAX™ qPCR and culture on at least a single sample, were PCR and culture negative. Tissue samples collected from rhinoceros number 10 were PCR negative, but culture positive for *M. bovis* (Supplementary Table 2). All the remaining culture-negative tissue and BAL samples tested were PCR negative (Supplementary Table 2).

Of the 31 antemortem and postmortem samples (6 tissues, 13 BALF and 12 TW) collected, from 13 African elephants, *M. tuberculosis* and *M. bovis* infection was confirmed by mycobacterial culture in 2 elephants; a zoo elephant and 1 free-ranging elephant from a bTB endemic area, respectively (Table 1 and Supplementary Table 2). The VetMAX™ qPCR confirmed infection in all culture-positive samples from these animals. Antemortem samples collected from the *M. tuberculosis*-infected elephant (two BALF and one TW) were all PCR positive, as well as a postmortem lung tissue sample. Only postmortem samples were available from the *M. bovis*-infected elephant. Five lymph node samples were all PCR positive, although only 1 lymph node was culture positive. The remaining 22 antemortem BALF and TW samples collected from 11 apparently healthy immobilized free-ranging elephants were all PCR and culture negative (Table 1 and Supplementary Table 2).

**Table 1** VetMAX™ PCR Kit test results for tissue and respiratory samples collected from African buffaloes (n = 48), white rhinoceros (n=10) and African elephants (n=13) with known mycobacterial culture results.

Species	Sample type	Total samples tested (collected from)	Confirmed MTBC infected animals <sup>a</sup> (n = 26 buffaloes, n = 6 rhinos & n = 2 elephants)				Confirmed MTBC uninfected animals (n = 22 buffaloes, n = 4 rhinos & n = 11 elephants)	
			Culture positive samples		Culture negative samples		Culture negative samples	
			PCR Positive <sup>b</sup>	PCR Negative	PCR Positive	PCR Negative	PCR Positive	PCR Negative
African buffaloes <sup>c</sup>	Lymph nodes(LNs)	n = 47 (LN pools from 26 <i>M. bovis</i> -positive- and 22 uninfected animals)	25	0	0	0	12	10
	Lung	n = 1 (single <i>M. bovis</i> -positive animal)	1	0	0	0	0	0
White rhinoceros <sup>d</sup>	LNs	n = 23 (all <i>M. bovis</i> -positive animals)	9	11	2	1	0	0
	Lung	n = 6 (each <i>M. bovis</i> -positive animal)	3	3	0	0	0	0
	BALF	n = 5 (all <i>M. bovis</i> -uninfected animals)	0	0	0	0	0	5
African elephants <sup>e</sup>	LNs	n = 4 (single <i>M. bovis</i> -positive animal)	1	0	3	0	0	0
	Lung	n = 2 (single <i>M. bovis</i> - and single <i>M. tuberculosis</i> positive animal)	1	0	1	0	0	0
	BALF	n = 13 (single <i>M. tuberculosis</i> -positive and all uninfected animals)	2	0	0	0	0	11
	TW	n = 12 (single <i>M. tuberculosis</i> -positive and all uninfected animals)	1	0	0	0	0	11

<sup>a</sup>Confirmed MTB infected: *M. bovis*- and *M. tuberculosis*-positive by mycobacterial culture and strain typing (22).

<sup>b</sup>VetMAX PCR Kit is a real-time PCR assay that targets the IS6110 insertion element sequence present in the genome of *M. tuberculosis*, *M. bovis*, *M. africanum*, *M. Microti*, *M. canetti*, *M. caprea* and *M. pinnipedii*.

<sup>c</sup>All test-positive African buffaloes were euthanized on either or the single intradermal comparative tuberculin test result and an interferon-gamma (IFN- $\gamma$ ) release assay result.

<sup>d</sup>African rhinoceros BALF samples tested were collected from four individual animals suspected of *M. bovis* infection, but not confirmed.

<sup>e</sup>African elephant PCR positive/ culture positive BALFs, TW and lung sample were collected from one elephant confirmed to be infected with *Mycobacterium tuberculosis*, one PCR positive/ culture positive LN samples were collected from one *M. bovis* confirmed elephant.

## Discussion

This study is the first to report the use of the VetMAX<sup>TM</sup> MTBC qPCR Kit, a nucleic acid amplification test, for detection of MTBC infection in African wildlife species.

When the VetMAX<sup>TM</sup> MTBC qPCR kit were applied to a bTB endemic African buffalo population, animals positive for the PCR but with negative culture results were likely to be truly infected (Table 1). This observation was probably due the imperfect sensitivity of mycobacterial culture (11). Moreover, for skin and/or cytokine assay positive animals, our findings suggest that PCR may be a better gold standard for determining MTBC infection status when assessing other diagnostic test performance (14).

Furthermore, for all known MTBC-infected rhinoceros, elephants and buffaloes, the observed 18 PCR-positive, mycobacterial culture-negative results (Table 1), may be due to the increased sensitivity of PCR along with the sub-optimal performance of

mycobacterial culture (11,12). These results may also be explained by the presence of residual DNA from dead mycobacteria in these samples, which would lead to a false positive PCR result (16). However, the culture and PCR negative results in the presumed uninfected rhinoceros and elephants also suggests that the VetMAX™ MTBC qPCR Kit has good specificity, although this will need to be investigated in future studies.

The discordant results, due to failure of the VetMAX™ MTBC PCR kit to detect DNA in *M. bovis* culture-positive tissue samples, may be due to DNA degradation, as a result of several freeze-thaw cycles of the stored homogenates (17). The DNA was extracted from stored tissues that had been previously frozen and thawed, whereas the samples were processed for culture after only one freeze-thaw cycle. However, DNA degradation may have been avoided if homogenized tissue samples, destined for DNA extractions, were stabilized with an appropriate media prior to freezing, or the best option would be to perform DNA extraction on fresh samples.

Based on preliminary data showing that this qPCR can detect MTBC DNA across a range of wildlife species and sample types, this kit appears to be a promising new tool for postmortem TB diagnosis as well as screening antemortem clinical samples. Since this kit can be used across species, it is especially important for wildlife where only limited host-based diagnostic tests are available. It is important to note that the technique appears to be sensitive, although the limitations of this pilot study included limited sample size and freeze-thawing of samples prior to qPCR testing.

## Conclusions

Our findings support the use of the VetMAX™ MTBC PCR Kit for the successful amplification of the insertion element IS6110 from mycobacterial DNA isolated from tissue and respiratory tract samples obtained from *M. bovis*-infected African buffaloes, *M. bovis*-infected white rhinoceros, and *M. tuberculosis*- and *M. bovis*-infected African elephants. Furthermore, this qPCR could provide a new tool for the detection of TB in various sample types from African wildlife species.

## Methods

### Animals

Archived tissue samples from forty-eight African buffaloes from bTB endemic areas, tested between 2016 and 2018 during a bTB test-and-cull programme, were selected for this study. All buffaloes were tested using the single comparative intradermal tuberculin test (SCITT), two separate interferon-gamma (IFN- $\gamma$ ) release assays (IGRAs) and a QuantiFERON TB Gold In-tube (QFT) IFN- $\gamma$  induced protein of 10 kDa (IP-10) release assay (IPRA), as previously described (18,19). All test-positive (i.e., positive on any test) buffaloes were shot by their respective park managers and transported to the abattoir for further processing. Thereafter, tissue samples were opportunistically collected at necropsy for mycobacterial culture and subsequent DNA extractions as in the case of this study. Number of animals selected for bTB testing were done by park management depending on funding available and not by the authors of this study.

Similarly, in 2018 in KNP, antemortem BALF samples were opportunistically collected from four boma-confined white rhinoceros and postmortem tissue samples from six poached rhinoceros. Trunk wash (TW) and BALF samples were also opportunistically collected from 12 chemically immobilized African elephants that included 11 free-ranging animals from KNP and one zoo elephant (South Africa) suspected of being MTBC-infected. In addition, tissue samples were collected from the same zoo elephant after natural death and opportunistically from an elephant carcass discovered in Kruger National Park. All rhinoceros and elephant tissue and respiratory tract samples were processed for DNA extractions and mycobacterial culture. Rhinoceros and elephant immobilizations for BALF and TWs collection were conducted by veterinarians from KNP Veterinary Wildlife Services (South African National Parks) according to their protocols, animals were carefully monitored during this process and all successfully revived. All animals reported in this study were opportunistically sampled during routine immobilizations and not specifically for this study.

### Sample collection and mycobacterial culture

All animal carcasses included in this study were examined for gross TB lesions. Tissue samples with TB-consistent lesions were collected and stored frozen separately. For animals with no visible lesions, samples of head (mandibular, retropharyngeal) and thoracic (mediastinal and tracheobronchial) lymph nodes were pooled by anatomical site. Endoscopic bronchoalveolar lavages were performed on chemically immobilized rhinoceros and elephants as previously described to collect antemortem samples (20,21). Additionally, TW samples were obtained from immobilized elephants by flushing 250 ml of sterile saline into each nostril, separately, elevating and lowering the trunk for approximately 30-40 s and aspirating the fluid into a sterile 500 ml collection chamber. Respiratory samples collected in 50 ml sterile falcon tubes were concentrated by centrifugation at 2000 x g for 30 min and decanting 46 ml of the supernatant. The remaining 4 ml with resuspended pellet were stored at - 20°C. All tissue and respiratory samples were processed for mycobacterial culture using Mycobacterial Growth Indicator Tubes (MGIT™) and the BACTEC™ MGIT™ 960 Mycobacterial Detection System (Becton Dickinson, Franklin Lakes, NJ, USA), as previously described (18). All Ziehl-Neelsen positively stained bacterial cultures were genetically speciated by PCR, as previously described (22).

### VetMAX™ MTBC PCR Kit

After mechanical lysis of all tissue, BALF, and TW samples, DNA was extracted from each sample using the QIAamp® DNA mini kit (Qiagen, Venlo, Limburg, Netherlands) as prescribed by the commercial VetMAX™ MTBC qPCR kit. Extracted DNA (5 µl) was mixed with 20 µl of reaction mix from the PCR kit and the reactions carried out in duplicate using the CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad Laboratories, Hercules, CA, USA). The PCR reactions were done at 50°C for 2 min (1 cycle), followed by one cycle at 95°C for 10 min, 40 cycles of 15 s at 95°C, and 1 min at 60°C. The experiment was conducted once due to reagent availability. Results were interpreted as negative ( $C_t > 40$ ), positive ( $C_t < 40$ ) or invalid, following the manufacturer's recommendations and by comparison with negative and positive qPCR controls. All samples were de-identified prior to DNA extractions to ensure unbiased interpretation of qPCR results and tested in batches of five to avoid possible cross contaminations.

### Statistical analysis

All positive and negative VetMAX™ qPCR and mycobacterial culture results were reported as proportions of the total number of samples tested. All results for the VetMAX™ qPCR, average  $C_t$  values and mycobacterial culture results for each sample tested from individual animals are provided in Supplementary Tables 1 and 2.

## References

1. Müller B, Dürr S, Alonso S, Hattendorf J, Laise CJM, Parsons SDC, et al. Zoonotic *Mycobacterium bovis*-induced Tuberculosis in humans - Volume 19, Number 6—June 2013 - Emerging Infectious Diseases journal - CDC. 2019 Aug 15 [cited 2019 Aug 15]; Available from: [https://wwwncdc.cdc.gov/eid/article/19/6/12-0543\\_article](https://wwwncdc.cdc.gov/eid/article/19/6/12-0543_article)
2. Dippenaar A, Parsons SDC, Miller MA, Hlokwe T, Gey van Pittius NC, Adroub SA, et al. Progenitor strain introduction of *Mycobacterium bovis* at the wildlife-livestock interface can lead to clonal expansion of the disease in a single ecosystem. *Infect Genet Evol J Mol Epidemiol Evol Genet Infect Dis*. 2017 Jul; 51:235–8.
3. Hlokwe TM, van Helden P, Michel AL. Evidence of increasing intra and inter-species transmission of *Mycobacterium bovis* in South Africa: are we losing the battle? *Prev Vet Med*. 2014 Jul 1;115(1–2):10–7.
4. Otchere ID, Tonder AJ van, Asante-Poku A, Sánchez-Busó L, Coscollá M, Osei-Wusu S, et al. Molecular epidemiology and whole genome sequencing analysis of clinical *Mycobacterium bovis* from Ghana. *bioRxiv*. 2018 Dec 6;489203.
5. Parsons SDC, Menezes AM, Cooper D, Walzl G, Warren RM, van Helden PD. Development of a diagnostic gene expression assay for tuberculosis and its use under field conditions in African buffaloes (*Syncerus caffer*). *Vet Immunol Immunopathol*. 2012 Aug 15;148(3–4):337–42.
6. Hlokwe TM, Said H, Gcebe N. *Mycobacterium tuberculosis* infection in cattle from the Eastern Cape Province of South Africa. *BMC Vet Res*. 2017 Oct 10;13. Available from: <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC5635548/>
7. Miller M, Michel A, van Helden P, Buss P. Tuberculosis in rhinoceros: An underrecognized threat? *Transbound Emerg Dis*. 2017 Aug;64(4):1071–8.

8. Miller MA, Buss P, Roos EO, Hausler G, Dippenaar A, Mitchell E, et al. Fatal Tuberculosis in a free-ranging African elephant and One Health implications of human pathogens in wildlife. *Front Vet Sci* [Internet]. 2019 [cited 2019 Feb 8];6. Available from: <https://www.frontiersin.org/articles/10.3389/fvets.2019.00018/full>
9. Miller MA, Buss P, Parsons SDC, Roos E, Chileshe J, Goosen WJ, et al. Conservation of White rhinoceroses threatened by Bovine tuberculosis, South Africa, 2016–2017. *Emerg Infect Dis*. 2018 Dec;24(12):2373–5.
10. Chileshe J, Goosen WJ, Buss PE, van Helden PD, Warren R, Parsons SDC, et al. A commercial ELISA for detection of interferon-gamma in white rhinoceros. *J Vet Diagn Investig Off Publ Am Assoc Vet Lab Diagn Inc*. 2019 Apr 11;1040638719843955.
11. de la Rua-Domenech R, Goodchild AT, Vordermeier HM, Hewinson RG, Christiansen KH, Clifton-Hadley RS. Ante mortem diagnosis of tuberculosis in cattle: a review of the tuberculin tests, gamma-interferon assay and other ancillary diagnostic techniques. *Res Vet Sci*. 2006 Oct;81(2):190–210.
12. Lorente-Leal V, Liandris E, Castellanos E, Bezos J, Domínguez L, de Juan L, et al. Validation of a real-time PCR for the detection of *Mycobacterium tuberculosis* complex members in bovine tissue samples. *Front Vet Sci* [Internet]. 2019 [cited 2019 Apr 2];6. Available from: <https://www.frontiersin.org/articles/10.3389/fvets.2019.00061/full>
13. Costa P, Botelho A, Couto I, Viveiros M, Inácio J. Standing of nucleic acid testing strategies in veterinary diagnosis laboratories to uncover *Mycobacterium tuberculosis* complex members. *Front Mol Biosci*. 2014;1:16.
14. Chakravorty S, Simmons AM, Rowneki M, Parmar H, Cao Y, Ryan J, et al. The new Xpert MTB/RIF Ultra: Improving detection of *Mycobacterium tuberculosis* and resistance to Rifampin in an assay suitable for point-of-care testing. *mBio*. 2017 Sep 6;8(4):e00812-17.
15. Dorman SE, Schumacher SG, Alland D, Nabeta P, Armstrong DT, King B, et al. Xpert MTB/RIF Ultra for detection of *Mycobacterium tuberculosis* and rifampicin resistance: a prospective multicentre diagnostic accuracy study. *Lancet Infect Dis*. 2018 Jan 1;18(1):76–84.
16. Theron G, Zijenah L, Chanda D, Clowes P, Rachow A, Lesosky M, et al. Feasibility, accuracy, and clinical effect of point-of-care Xpert MTB/RIF testing for tuberculosis in primary-care settings in Africa: a multicentre, randomised, controlled trial. *Lancet Lond Engl*. 2014 Feb 1;383(9915):424–35.
17. Shao W, Khin S, Kopp WC. Characterization of effect of repeated freeze and thaw cycles on stability of genomic DNA using pulsed field gel electrophoresis. *Biopreservation Biobanking*. 2012 Feb;10(1):4–11.
18. Goosen WJ, Miller MA, Chegou NN, Cooper D, Warren RM, van Helden PD, et al. Agreement between assays of cell-mediated immunity utilizing *Mycobacterium bovis*-specific antigens for the diagnosis of tuberculosis in African buffaloes (*Syncerus caffer*). *Vet Immunol Immunopathol*. 2014 Jul 15;160(1–2):133–8.
19. Goosen WJ, Cooper D, Miller MA, van Helden PD, Parsons SDC. IP-10 Is a sensitive biomarker of antigen recognition in whole-blood stimulation assays used for the diagnosis of *Mycobacterium bovis* Infection in African buffaloes (*Syncerus caffer*). *Clin Vaccine Immunol CVI*. 2015 Aug;22(8):974–8.
20. Michel AL, Lane EP, de Klerk-Lorist L-M, Hofmeyr M, van der Heijden EMDL, Botha L, et al. Experimental *Mycobacterium bovis* infection in three white rhinoceroses (*Ceratotherium simum*): Susceptibility, clinical and anatomical pathology. *PloS One*. 2017;12(7):e0179943.
21. Parsons SDC, Morar-Leather D, Buss P, Hofmeyr J, McFadyen R, Rutten VPMG, et al. The kinetics of the humoral and Interferon-gamma immune responses to experimental *Mycobacterium bovis* infection in the White rhinoceros (*Ceratotherium simum*). *Front Immunol* [Internet]. 2017 Dec 22 [cited 2019 Mar 13];8. Available from: <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC5743672/>
22. Warren RM, Gey van Pittius NC, Barnard M, Hesseling A, Engelke E, de Kock M, et al. Differentiation of *Mycobacterium tuberculosis* complex by PCR amplification of genomic regions of difference. *Int J Tuberc Lung Dis Off J Int Union Tuberc Lung Dis*. 2006 Jul;10(7):818–22.

## Abbreviations

1. MTBC: *Mycobacterium tuberculosis* complex
2. *bovis*: *Mycobacterium bovis*

3. *tuberculosis*: *Mycobacterium tuberculosis*
4. bTB: bovine tuberculosis
5. TB: tuberculosis
6. KNP: Kruger National Park
7. CMI: cell-mediated immunity
8. BALF: bronchoalveolar lavage fluid
9. TW: trunk wash
10. RT-qPCR: real time quantitative polymerase chain reaction
11. IS6110: insertion element 6110
12. IFN- $\gamma$ : interferon gamma
13. IP-10: interferon gamma-induced protein of 10kDa
14. SCITT: single comparative intradermal tuberculin test
15. IGRA: IFN- $\gamma$  release assay
16. IPRA: IP-10 release assay
17. QFT: QuantiFERON<sup>®</sup> TB Gold In-tube assay
18. MGIT: Mycobacterial growth indicator tuber
19. ZN stains: Ziehl-Neelsen staining

## Declarations

### Ethics approval

Ethical approval for this project was granted by Stellenbosch University Animal Care and Use Committee (ACU-2018-0966; ACU-2019-6308 and ACUACU-2019-9081) and section 20 research permits issued by the Department of Agriculture, Forestry and Fisheries (Ref: 12/11/1/7/2; 12/11/1/7/5 and 12/11/1/7/6).

### Consent for publication

Not applicable.

### Availability of data and materials

The raw data generated during the current study are available from the corresponding author on reasonable request. Moreover, all results for each animal are available in the accompanying supplementary tables.

### Competing interest

The authors declare that they have no competing interests.

### Funding

This work was only supported financially by 1) the Harry Crossley Foundation, 2) South African government through the South African Medical Research Council and the National Research Foundation South African Research Chair Initiative (grant #86949), 3) ThermoFisher, and 4) American Association of Zoo Veterinarians Wild Animal Health Fund (S005651). The study design, content and conclusions are the sole responsibility of the authors and does not represent the official views of the funders. The funder's only role was for financial assistance.

### Author's contributions

WJG and MM designed the study. WJG performed the experiments. WJG, MM and LK analysed the data. WJG, MM, LK, RMW, Pvh, BS and SDCP wrote the manuscript. PB and DC assisted with animal sampling. All authors read and approved the final

manuscript.

## Acknowledgements

We thank Alicia and Warren McCall, Dr. Rowan Leeming, JP van Heerden, the Game Capture staff from Kruger National Park and Enzemvelo Wildlife, Guy Hausler, Leana Rossouw, Tebogo Manamela, and Drs. Jennie Hewlett and Robert Campbell from National Zoological Gardens of South Africa.

## Author's information

<sup>a</sup>DST-NRF Centre of Excellence for Biomedical Tuberculosis Research; South African Medical Research Council Centre for Tuberculosis Research; Division of Molecular Biology and Human Genetics, Faculty of Medicine and Health Sciences, Stellenbosch University, Cape Town, South Africa.

<sup>b</sup>Veterinary Wildlife Services, South African National Parks, Kruger National Park

<sup>c</sup>Enzemvelo KZN Wildlife, PO Box 25, Mtubatuba 3935, South Africa

<sup>d</sup>Thermo Fisher Scientific; Prionics AG; Wagistrasse 27A; Schlieren-Zurich; Switzerland

## Supplementary Files

This is a list of supplementary files associated with this preprint. Click to download.

- [Table1VetMAXTMPCKittestresults..docx](#)
- [VetMAXMTBCPCRKitWJG2020BMCV2.docx](#)