

# Wideband MR Echo Planar Imaging: A High Temporal/Spatial Resolution Neuro-architecture Study of Rodent Brain

Po-Wei Cheng

Graduate Institute of Biomedical Electronic and Bioinformatics, National Taiwan University

Tzi-Dar Chiueh

Department of Electrical Engineering, National Taiwan University

Jyh-Hong Chen (✉ [jhchen@ntu.edu.tw](mailto:jhchen@ntu.edu.tw))

Department of Electrical Engineering, National Taiwan University

---

## Research Article

**Keywords:** Latest simultaneous multi-slice (SMS) methods, parallel imaging technique, precision neuro-architecture studies

**Posted Date:** June 1st, 2021

**DOI:** <https://doi.org/10.21203/rs.3.rs-413062/v2>

**License:**  This work is licensed under a Creative Commons Attribution 4.0 International License.

[Read Full License](#)

---

# Wideband MR Echo Planar Imaging: A High Temporal/Spatial Resolution Neuro-architecture study of Rodent Brain

Po-Wei Cheng<sup>1</sup>, Tzi-Dar Chiueh<sup>2</sup>, Jyh-Horng Chen<sup>1,2,3\*</sup>

<sup>1</sup> Graduate Institute of Biomedical Electronic and Bioinformatics, National Taiwan University, Taipei, Taiwan

<sup>2</sup> Department of Electrical Engineering, National Taiwan University, Taipei, Taiwan

<sup>3</sup> Interdisciplinary MRI/MRS Lab, Department of Electrical Engineering, National Taiwan University, Taipei, Taiwan

\* Corresponding authors: Author Three (jhchen@ntu.edu.tw).

## Abstract

**Latest simultaneous multi-slice (SMS) methods greatly benefit MR efficiency for recent studies using parallel imaging technique. However, these methods are limited by the requirement of array coils. The proposed Coherent Wideband method, which employs an extended field of view to separate multiple excited slices, can be applied to any existing MRI instrument, even those without array coils. In this study, the Coherent Wideband echo-planar imaging method was implemented on 7T animal MRI to exhibit comprehensive enhancements in neuro-architecture, including diffusion tensor imaging (DTI) and functional MR studies (fMRI). Under the same scan time, the time-saving effect can be manipulated to increase the number of averages for DTI SNR improvement, reducing fractional anisotropy difference by 56.9% (from 0.072 to 0.041) and the deviation angle by 64% (from 25.3° to 16.2°). In summary, Coherent Wideband EPI will provide faster, higher resolution, thinner slice, or higher SNR imaging for precision neuro-architecture studies.**

## Introduction

Magnetic resonance imaging (MRI)<sup>1</sup> technology is a non-invasive diagnostic tool that reveals not only changes in cerebral blood oxygen saturation<sup>2,3</sup> but also pathways in nerve fibers.<sup>4-6</sup> However, compared with computed tomography and positron emission tomography, MRI has a markedly lower time resolution.

Consequently, numerous acceleration techniques were born in pursuit of greater MRI throughput. Simultaneous multi-slice (SMS) imaging is an acceleration technique which enhances imaging efficiency by simultaneously exciting and acquiring multiple slices. Müller first proposed the concept in 1988 and used the Fourier shift theorem to develop multi-frequency selective radiofrequency pulses for SMS<sup>7</sup>. This led to the proposal of Hadamard-encoded RF pulses, which Souza et al. utilized to separate simultaneously acquired slices.<sup>8</sup> Glover et al. took a similar approach, utilizing encoded RF pulses in their phase-offset multi-planar (POMP) technique.<sup>9</sup> These earliest SMS methods didn't shorten MR scan time due to the extra phase-encoding steps or excitations needed in their design. However, they facilitated the development of future SMS techniques in the next few decades.

Modern SMS methods can be divided into two categories according to their hardware dependency. The majority of SMS methods belong to the first category, which utilizes additional hardware to provide

necessary spatial information to separate simultaneously acquired slices. Multi-channel coils were first proposed and used by Larkman to acquire multiple excitations.<sup>10</sup> Afterward, Breuer used a combination of RF and multi-channel coils for multi-slice imaging.<sup>11</sup> Furthermore, Setsompop et al. combined controlled aliasing in parallel imaging with echo-planar imaging to match multi-excited image signals with standard signals.<sup>12,13</sup> These SMS techniques hugely advanced the research in fields of functional MRI (fMRI) and diffusion tensor imaging (DTI)<sup>14,15</sup>. Numerous studies successfully proved the benefits of applying SMS techniques instead of other techniques, such as increasing analysis bandwidth<sup>16</sup> or enhances image resolution.<sup>17</sup>

However, all aforementioned SMS methods suffer from hardware-induced artifact problems like g-factors and slice leakage due to the usage of multi-channel coils.<sup>18</sup> Since small animal pre-clinical apparatuses generally possess fewer coils, the noise will be even more severe. From this aspect, non-hardware based SMS approaches are essential to provide neural dynamics in pre-clinical studies.

Weaver *et al.* first proposed the approach to achieve simultaneous multi-slice imaging acceleration without extra hardware in 1988.<sup>19</sup> The multi-excited signals can be separated by applying an extra gradient during spatial encoding, but this addition results in severe image blur. Therefore, Wu *et al.* developed the Multi-frequency excited Wideband (ME-Wideband)<sup>20</sup> method on gradient sequences for blur mitigation, improving Wideband image quality to the standard level. Since most current SMS techniques use additional hardware, non-hardware based SMS techniques are more commonly referred to as “Wideband” techniques to emphasize the difference.<sup>21,22</sup>

For dynamic functional studies, Setsompop et al. applied Weaver’s technique to EPI and referred to it as Blipped-Wideband.<sup>12</sup> Naturally, Blipped-wideband suffers from the same blurring problem. For blur mitigation, the ME-wideband technique is applied to EPI in this paper. We adapt ME-wideband to EPI sequence, the “Coherent Wideband” technique is proposed to optimize wideband EPI images. Coherent Wideband, consisting of novel refocusing gradient sequence and precise phase alignment, provides faster, higher resolution, thinner slice imaging, or higher signal-to-noise pre-clinical imaging for precision neuro-architecture studies. Experimental results prove the advantages of Coherent Wideband against previous Wideband methods and demonstrate the benefits of this technique in fields of pre-clinical DTI and fMRI researches.

## Materials and Methods

### Pulse sequence and imaging parameters

This paper proposes the “Coherent Wideband” technique, consisting of a novel refocusing gradient sequence and precise phase alignment. This study introduced a method that employs expanded field of view (FOV) to acquire image signals from simultaneously-excited slices without extra hardware, then utilizes separation gradients so that image signals undergo phase shifts. Although additionally applied separation gradient and expanded FOV can separate image signal from different slices on the image domain, they also cause residual artifacts or “N/2 ghosting” artifacts. To eliminate these artifacts, an accurate phase alignment method is introduced by this study.

### Coherent Wideband EPI Sequence

For EPI sampling, extra separation gradient pulses were applied along the z-direction to achieve phase shift in signals and corresponding image separation. For two-fold excitation (acceleration factor,  $W=2$ ), each  $k$ -line exhibits an extra phase of  $\pi/2$ , which results in a FOV/2 shift on the image.


According to the ME-wideband equation by Wu, signal strength and phase vary by this separation gradient applied as describe in Eq1&2.

$$S'(k_x, k_y) = S(k_x, k_y) \times \left( \int_{z_1}^{z_2} \exp(i\gamma G_z \tau) dz \right) \quad \text{Eq.1}$$

$$\left( \int_{z_1}^{z_2} \exp(i\gamma G_z \tau) dz \right) = (z_2 - z_1) \times \text{sinc}\left(\gamma G_z \tau \left(\frac{z_2 - z_1}{2}\right)\right) \times \exp\left(i\gamma G_z \tau \left(\frac{z_2 + z_1}{2}\right)\right) \quad \text{Eq.2}$$

A gradient is required to separate two adjacent planes, which further reduces the signal strength in wideband MRI. Blipped-CAIPI method<sup>23</sup> and ME-Wideband<sup>20</sup> have suggested that refocusing gradients can be applied to recover signal strength.

Fig.1 shows multiple wideband EPI sequence designs, including blipped wideband EPI, ME-wideband EPI, and Coherent Wideband EPI. ME-wideband EPI is the EPI version of ME-Wideband, while Coherent Wideband EPI is its optimization. In blipped wideband EPI, accumulated separation gradient causes severe signal attenuation. While ME-wideband EPI partially mitigates such attenuation, the rapidly switching refocusing gradient leads to additional eddy current artifact. The optimized Coherent Wideband design keeps the isochromatic spin in phase in each pair of separation gradients, minimizing signal attenuation, phase deviation, and eddy current artifacts. (See Fig.1 (d), (e) & (f)) The phase plots of Coherent Wideband shows the phase shifts between 0 and  $\pi/2$  which causes an N/2 shift on the image domain. The following section will demonstrate how an accurate phase alignment process addresses this issue.




**Fig. 1. Pulse sequences of multiple wideband EPI**

(a) Blipped Wideband EPI Sequence with constant separation gradients. (b) ME Wideband EPI Sequence with refocusing gradient. (c) Coherent Wideband EPI sequence with bipolar refocusing gradient. (d) Phase plots of multiple EPI sequences. (e) Signal Magnitude of various EPI sequences. (f) Zoom in View of Signal Magnitude

### Zero phase $N/2$ artifact mitigation: EPI Reconstruction

EPI images suffer from  $N/2$  artifacts due to rapidly switching gradient and resulting odd-even echo phase errors. Since these artifacts overlaps with multi-slice images in Wideband EPI applications,  $N/2$  ghost removal is the priority of our image reconstruction. Proper phase corrections can be done to reduce  $N/2$  ghost by measuring the gradient data set in Fig. 2a & 2b. <sup>24</sup> However, ghost artifacts are still visible due to k-space signal asymmetric. The zero-phase compensation method based on Ordidge *et al.* <sup>25</sup> improves the phase adjustment process, facilitating the symmetry of k-space signal and reducing the ghost factor from 6.3% to 3.1% (Fig. 2c). Nevertheless, for the case of wideband EPI, the phase difference caused by additional separation gradient should also be considered. The next section will elaborate on this extra phase alignment process.




**Fig. 2. EPI Phase Alignment.**

Multiple  $N/2$  ghost correction method is implemented in (a), (b) & (c). The uncorrected  $N/2$  ghost artifact is shown in (a.3). The artifact is partially removed by the linear correction method in (b.3), and thoroughly removed by the zero-phase compensation method in (c.3).

### Slice Dependent Phase Alignment: Wideband Reconstruction

The Wideband technique utilizes additional separation gradients to separate simultaneously-excited images. These gradients introduce different amount of phase to signals depending on the spatial location of the excited slices, which also contributes to N/2 artifacts as shown in Fig.3 (a). Therefore, this extra phase must be measured and compensated to ensure image quality. Ideally, this phase (Fig.3 (b)) can be calculated with the given sequence design, and phase correction can be done accordingly to remove ghost artifacts. However, the actual phase value is affected by the inhomogeneity of magnetic field and eddy currents, so remnant artifacts could still be observed after phase correction with theoretical phase values, as shown in Fig. 3(c). To ensure precise phase correction for Wideband image reconstruction, the gradient measurement method proposed by Beaumont<sup>26</sup> is applied in this study to determine the magnitude and phase of signals, as shown in Fig. 3(d). As Fig. 3(e) demonstrates, the ghost artifacts are no longer visible after phase correction with measured phase values.



**Fig. 3. Slice Dependent Phase Alignment.**

Separation gradients induce phase shift in k-space which leads to N/2 artifacts on the images shown in (a). The artifacts were partially corrected by theoretical phase as shown in (b) & (c). By applying the measured phase, precise offset compensation for Wideband images was ensured as shown in (d) & (e).

### Animals

In this study, pre-clinical scans were performed on the male healthy Sprague Dawley rat model (200-300g, N=6 in this study) from National Taiwan University Laboratory Animal Center. Pre-clinical experimental procedures were approved by the Institutional Animal Care and Use Committee of National Taiwan University. All animal experimentation was done with institutional approval and followed NTU guidelines. This study is conducted in accordance to ARRIVE guidelines.

### MR Parameters

All experimental data used in this study was obtained on the Bruker 7T MRI system by a volume coil (Biospec 70/30, PV5.1, Bruker BioSpin MRI GmbH, Ettlingen, Germany). All experimental protocols were approved by Use Committee of National Taiwan University.

- *Rodent brain structure study: High temporal resolution*

Coherence Wideband EPI achieves 2-fold acceleration for rodent brain scan. The scan parameters are as follows: For conventional Spin-Echo EPI, TR/TE = 3000/75 ms, matrix size = 96x192, slice thickness = 1mm, resolution = 0.26 x 0.26 mm<sup>2</sup>, bandwidth = 250 kHz, the number of averages = 120, and the total time is 6 mins. For 2-fold acceleration Coherent Wideband EPI, TR / TE = 1500 /75 ms, cutting the scan time to 3 min.

- *Rodent brain structure study: High spatial resolution*

Another usage of Coherent Wideband EPI is enhancing image spatial resolution under the same scan duration. The scan parameters are as follows: For traditional Spin-Echo EPI, TR/TE = 1000/75 ms, matrix size = **80x80**, slice thickness = 1mm, resolution = 0.31 x 0.31 mm<sup>2</sup>, bandwidth = 500 kHz, the number of averages = 40, and the total time is 6 m 40 s. For 2X-spatial resolution Coherent Wideband EPI, matrix size = **160x160**, TR/TE = 500/75, and the resulting scan duration is also 6 min 40 s.

- *Rodent diffusion Tensor image: High temporal resolution*

The parameters for EPI-based Spin-Echo DTI are as follows: FOV = 2.5x5.0 cm<sup>2</sup>, slice thickness = 1mm, matrix Size = 96x192, TE =75ms, and the number of averages = 40. TR = 5000 ms for standard EPI and 2500 ms for Coherent Wideband EPI, so the respective scan time is 43 min 20 s and 21 min 40 s.

- *Rodent diffusion Tensor image: SNR enhancement*

The time-saving advantage of the Coherent Wideband technique can also elevate image SNR by doubling the number of average for precision medicine study. In this study, the Spin-Echo DTI scan consists of 1 null-image and 12-direction diffusion gradients. The scan parameters are as follows: FOV=2.5x5.0cm<sup>2</sup>, slice thickness = 1mm, TE =75ms,  $\Delta/\delta = 8/3$ ms, b-value = 1500s/mm<sup>2</sup>. TR = 5000 ms for standard EPI and 2500 ms for Coherent Wideband EPI. Since Coherent Wideband EPI is two times faster, the number of average is increased from 40 to 80 **while total scan time remains constant at 43min 20s.**

- *Rodent diffusion Tensor image: Thinner Slice*

For higher Z-axis resolution, Coherent Wideband EPI provides thinner slice imaging to delineate finer details of brain neuro-architecture. The parameters are as follows: For conventional EPI, slice thickness = 1 mm, average number = 20, **number of slices = 12** and TR = 5000; for Coherent Wideband EPI, slice thickness = 0.5 mm, number of average = 20, **number of slices = 24** and TR = 2500. The duration of both scans is 32 m 30 s.

- *Rodent resting functional MRI: 1/2 short TR*

High sampling MR speed is one of the unmet demands for resting fMRI researches (N=6). For investigating the quick dynamic change in fMRI, wideband EPI is utilized to improve temporal resolution in each measurement. The TR of Coherent wideband is shortened from 2 s to 1 s, and the sampling point for each fMRI study is increased from **120 to 240** for better statistical significance of functional analysis under the same scan time. The alternative is to shorten the acquisition time for **120** points from 6min to 3min to fit the user's need.



## MR Analysis

### Diffusion Analysis

The k-space data of the wideband images were reconstructed in Matlab (Math Works, Natick, USA). The DSI Studio (<http://dsi-studio.labsolver.org>) was used to perform Diffusion Tensor analysis and calculating the eigenvalues for each pixel. After selecting different ROIs, their respective Fractional Anisotropy (FA), Mean Diffusivity (MD), and other values can be calculated to determine their directivity and diffusivity<sup>27,28</sup>, and then eigenvectors can be used to calculate the direction and angle of the pixel<sup>29</sup>.

### Reproducibility Study

Even with the same scan parameters and conditions, different repetitions of the same DTI scan may derive different FA and fiber orientation due to thermal and systematic noise<sup>29</sup>. Therefore, high SNR is required for DTI to reduce the variance of analysis results and consistently track neuro-structures. To observe the reproducibility rise for SNR-enhanced DTI results, these DTI scans will be repeated twice to calculate their FA difference and deviation angle between repetitions.

### Resting MRI analysis

Resting functional MR image analysis was performed using REST toolbox (<http://www.restfmri.net/forum/>), with SPM 12 (Wellcome Department of Cognitive Neurology, London, UK; <http://www.fil.ion.ucl.ac.uk/spm>). Default preprocessing pipeline was used to perform all the preprocessing steps, including slice-time correction, and then spatially smoothed. Physiological and other spurious noise were estimated and then removed together with the aforementioned motion artifacts by band-pass filtering (chosen band 0.001Hz-1Hz). The ROIs were defined and used as a SEED, and functional ROI connected to each Segment. For first-level analysis, the residual BOLD time course was extracted from each ROIs and Pearson's correlation coefficients were computed between each pair of ROIs.

## Data availability


The image data generated and presented in this study are available from the corresponding author on reasonable request.

## Results

The following figure presents phantom images acquired by standard spin-echo EPI and various Wideband-EPI ( $W = 2$ ) sequences in Fig.1. For image quality evaluation, two repeat measurements were performed for each set of images to calculate its SNR<sup>30</sup>, standard deviation (SD) value of different interest regions<sup>31</sup>, and SSIM<sup>32</sup>. The standard EPI phantom images were used as the benchmark.

Figure 4 indicates that the blipped wideband EPI, which uniformly exerts separation gradient, generated blurry images. The blur causes a drastic drop in SSIM, from 0.998 to about 0.58. The ME-Wideband EPI partially mitigated the blur improving SSIM to 0.93. However, the image results still show visible inhomogeneity as indicated by the high standard deviation value of different regions. By contrast, the Coherent EPI method, which applies bipolar separation gradients, resolved the inhomogeneity problem

as SD reduces from 138.1 to 113.1 and SSIM further improves to 0.995. In summary, among all wideband methods, Coherent Wideband EPI images are consistent with conventional results, as the SNR, SD, and SSIM values indicate.




**Fig. 4. Phantom results for multiple Wideband EPI methods.**

(a) Standard EPI results (TR/TE = 500/55 ms). Slices A & B (10 mm from center respectively) were separately excited and acquired. (b) Blipped Wideband EPI (TR/TE = 250/55ms) result. Slices A & B were simultaneously excited and acquired. (c) The ME-Wideband EPI partially mitigated the blur and showed visible inhomogeneity. (d) Coherent Wideband EPI results are consistent with conventional results.

### 2X Temporal Resolution Rodent Brain Results

Fig. 5 presents the 2X Coherent Wideband spin-echo EPI in-vivo rodent brain results acquired two times faster than a conventional EPI scan, while the image quality maintains the same as standard images. The slice position is coordinated with the rodent brain atlas.<sup>33</sup> Quantitative comparison by SNR and SSIM values also indicate similar image quality for both methods.



**Fig. 5. Conventional and 2-fold Temporal Wideband in vivo rodent brain study**


The In Vivo rodent brain images were acquired by a 7T animal MRI with (a) Conventional Spin-Echo EPI (b) Two-fold temporal resolution Coherent Spin-Echo EPI. Coherent Wideband EPI is 2 times faster while the image quality maintains the same as standard images.

### 2X Spatial Resolution Rodent Brain Result

Since the Coherent Wideband technique reduces the repetition of EPI scan, Coherent Wideband EPI only needs a TR of 500 ms to perform 10 scans while conventional EPI needs 1000 ms. As the result, the number of kx and ky can be doubled for Wideband EPI to achieve higher resolution under the same scan time.

The spatial resolution difference between conventional and wideband EPI scans can be compared. Under the same scan time, the matrix size is 80x80 points for standard EPI scan and 160x160 for Wideband (as shown in Fig.6).

Since the partial volume effect is reduced due to higher resolution, corresponding structural areas like blood vessels and cortical boundaries are more visible, and the demarcation for gray and white matter is clearer. Results show that the boundaries between these structures were vague due to the partial volume effect for low-resolution images. With the Wideband technique, the boundaries were much clearer due to spatial resolution enhancement.




**Fig. 6. Conventional and 2-fold spatial resolution Wideband rodent brain results.**

(a) For standard EPI, the image resolution was  $312 \times 312 \mu\text{m}^2$  and  $\text{TR}=1000\text{ms}$ , so the total scan time was 6 m 40 s. (b) 2X Coherent Wideband acquired double data points under the same scan time, improving the image resolution to  $156 \times 156 \mu\text{m}^2$ .

### **DTI Result**

Fig.7 presents high temporal resolution DTI results achieved by the Coherent Wideband spin-echo EPI technique. Under the same scan time, the time-saving effect can also be traded off against SNR improvement by increasing the number of averages, as demonstrated in Fig.8 to Fig.10. To reduce the partial volume effect of DTI tracking in fiber kissing and fiber crossing, the slices can also be thinned by half as shown in Fig.11 & Fig.12.



**Fig. 7. The 2-fold temporal resolution diffusion study.**


In the Coherent Wideband 2x temporal resolution diffusion experiment, the acquisition time of the sequence is reduced from 43 min 20 s to 21 min 40 s. The signal to noise and FA value of the 2X fast image kept well at the equal quality of the original one.

#### 2X-Average Number Coherent Wideband EPI DTI

In this experiment, a 1.4X SNR enhancement for 2X average Wideband images with the same acquisition time was observed as expected in Fig8. The comparison shows a slight FA value rise with a clear standard deviation decrease for 2-average Wideband enhanced DTI. To quantify the accuracy of diffusion MRI in axonal fiber orientation, each DTI scan was performed twice, and the fractional anisotropy (FA) difference and deviation angles (DA) for repeated scans were determined.

As Fig.8 demonstrates, 2X-Average Wideband DTI, resulting in higher SNR, the half width at half maximum (HWHM) of whole brain FA difference distribution dropped from 0.072 to 0.041, and the HWHM of fiber deviation angle histogram decreased from 25.3° to 16.2°. Both plots indicate that the reproducibility of FA and fiber angle improved along with the 1.4X SNR gain in Fig.9.


The conventional DTI tractography images and the SNR-enhanced Wideband DTI images are shown in Fig.10. With higher SNR, the Wideband DTI tractography shows a more robust tracking result in the external capsule (ec), corpus callosum (cc), and cingulum (cg) without noticeable erroneous forks.



**Fig. 8. Diffusion study of 2X-Average Coherent Wideband EPI images.**


In the 2X-Average Coherent Wideband MR diffusion experiment, the SNR was boosted 1.4 times from 76 to 105, as expected. The acquisition time of the sequence remains at 43min 20s.

Reproducibility study



**Fig. 9. Reproducibility analysis of FA and axonal fiber orientation**

The quantitative analysis of DTI result, including (a) fractional anisotropy (FA) difference and (b) deviation angle. Due to the 1.4x SNR gain, Wideband results have lower HWHM in both histograms, indicating better reproducibility. The half width at half maximum (HWHM) of whole brain FA difference distribution dropped from 0.072 to 0.041(56.9%), and the HWHM of fiber deviation angle histogram decreased from 25.3° to 16.2° (64%).



**Fig. 10. Diffusion Tracking of 2X-Average Coherent Wideband EPI.**


Images are magnified for the comparison between conventional EPI (SNR=76.1) and SNR enhanced Wideband EPI (SNR=105.2) tractography (2x scan numbers; same scan time). The Wideband result shows higher SNR and uniform tracking in the external capsule (ec), corpus callosum (cc), and cingulum (cg).

#### 2X Thinner slice Coherent Wideband EPI DTI

A low SNR diffusion tractography is susceptible to the partial volume effect which deteriorates analysis results. In this experiment, the SNR for DTI results was set to 55, a value lower than previous DTI studies, to demonstrate the partial volume effect for DTI analysis. Under these circumstances, a conventional DTI scan yields specious tracking results, particularly in the cc and cg regions. On contrary, the 2x thinner Coherent Wideband DTI tractography with 0.5mm thickness (Fig.11b) presents more detailed neuro-architecture and more accurate tracking results according to the atlas. As Fig.11 demonstrates, the demarcation of the axonal fiber orientation is also improved with thinner slices.


Fig. 12 is the coronal view of its DTI results, presenting the entire CC fiber tract. The FA maps illustrate the diffusivity, and the color maps represent the orientation of axons. As shown in the FA maps (Fig. 12a & 12b), the 2X-Thinner Wideband DTI result shows accurate fiber continuity, especially for the Hippocampal Commissure (HC) region. The consistency of fiber tracking for the CC region is also improved, as demonstrated in the color maps (Fig. 12c & 12d).





**Fig. 11. The diffusion tracking of 2X Thinner Wideband results.**

Conventional EPI DTI images (SNR=55) and thinner Wideband EPI DTI images (2x thinner slices; same scan time). The conventional low SNR DTI scan yields specious tracking results, particularly in the corpus callosum (cc), and cingulum (cg). The comparison shows that the Wideband result presents more details of neuro-architecture and consistent tracking in fimbria (fi), ec, cc, and cingulum (cg).




**Fig. 12. The 2x thinner wideband diffusion study.**

DTI FA maps (Fig. 12a & 12b) and DTI Color maps (Fig. 12c & 12d) are magnified to compare the conventional and thin-slice Wideband results. Since the partial volume effect is reduced, the precision of quantitative DTI analysis is the uniformity of tracking are enhanced in the Hippocampal Commissure (HC) and corpus callosum (cc).

## Resting Result


Since current resting fMRI analyses rely on regional locality, the accuracy of results highly depends on its temporal signal to noise ratio (tSNR)<sup>34</sup>, which means that more time points and higher temporal resolution can both improve the precision and stability of fMRI results. In the application results of resting fMRI, the Coherent Wideband gradient-echo EPI technique doubled the sampling rate to acquire more data points under the same scan time for better resting-fMRI analysis. The tSNR of Wideband resting-fMRI is enhanced in Cortex, CC, and Putamen (Fig.13). The tSNR was elevated by 44.2% from  $43.28 \pm 9.65$  to  $62.41 \pm 9.52$  by shortening the sampling interval and was boosted by 65.2% to  $70.9 \pm 9.22$  by doubling the scan number in the cortex area. This enhancement results directly from the reduction of thermal noise with high sampling point number and the reduction of physiological noise with high

sampling rate while using Coherent Wideband EPI. As the experimental result demonstrates, the Coherent Wideband EPI result shows a 38.3% (Fig.14a&14b) increase of associated activation region size in comparison with the traditional resting fMRI result due to the tSNR enhancement in the double-sampling case.



**Fig. 13. Temporal signal to noise ratio map.**

In the resting functional MRI experiment, the tSNR map was boosted by 65.2% to  $70.9 \pm 9.22$  by doubling the scan number, and was elevated from  $43.28 \pm 9.65$  to  $62.41 \pm 9.52$  (44.2%) by shorten the sampling interval in the cortex area. In the corpus callosum and Putamen, by doubling the scan number in the cortex area, tSNR increased by 46.5% and 42.1%; by shortening the sampling interval, tSNR increased from  $34.23 \pm 11.28$  and  $32.12 \pm 13.22$  to  $42.17 \pm 9.31$  (23.2%) and  $37.71 \pm 9.28$  (17.4%).




**Fig. 14. Resting-state functional connectivity.**



Resting-state functional connectivity with bandpass filtered at 0.01 and 1Hz. Group resting-state functional connectivity maps ( $p < 0.005$ ,  $N = 6$ ) using seed regions in the (a) right motor network, (b) right somatosensory network. The 2X Wideband EPI (Scan Number=240) shows increased (1.38 times as conventional scan) regional connectivity, but 2X Wideband EPI (Scan Number=120) reveals less (0.78 times as conventional scan) connections in the right motor network. In the right somatosensory network, 2X Wideband EPI (Scan Number=240) presents 1.34 times relative regions but 0.82 times regions are observed in 2X Wideband EPI (Scan Number=120).

## Discussion

This study presented multiple applications of the Coherent Wideband EPI technique and demonstrated its multifaceted benefits (Fig.15). For temporal resolution enhancement, Coherent Wideband EPI scans multiple slices simultaneously to reduce TR and thus total duration, providing results with double speed while the quality remains identical with the conventional results. This method can also obtain double data points under the same scan time to elevate the spatial resolution of rodent brain results. Regarding diffusion applications, under the same scan time, Coherent Wideband EPI can either decrease slice thickness to diminish the partial volume effect of DTI tractography or increase the number of averages to reduce FA difference and fiber deviation angle. As for resting fMRI, the proposed technique doubled the temporal resolution for a significant improvement in analysis results. Nonetheless, the Coherent Wideband EPI technique still has room for improvement in multiple aspects, which we'll respectively elaborate on in the following sections.



**Fig. 15. Multiple applications of the Wideband MRI.**

This study demonstrated multiple benefits of the Coherent Wideband EPI technique, including temporal/spatial resolution enhancement, SNR improvement, and accuracy study.

### TE limitations

To achieve a  $W$ -fold acceleration with Wideband, the imaging process acquires  $W$  times as many points as the conventional scan. These additional data points prolong the TE of Wideband scans and

exacerbate T2 decay, which could be more severe for the high field. To resolve this issue, the receiver bandwidth for Coherent Wideband EPI should be also increased to fully restore signal strength. Hence, the gradient specifications of MRI systems, including slew rate and gradient amplitude, need to be compatible with wider bandwidth for the highest possible W for Coherent Wideband. For the current Bruker 7T Biospec animal system, the maximum achievable W is 3. To achieve further acceleration rate solely with the Coherent Wideband technique, a gradient upgrade is required.

Nevertheless, the Coherent Wideband EPI technique has high compatibility with other acceleration methods, including Hadamard encoding and parallel imaging techniques. Therefore, further acceleration is achievable without further increasing TE or bandwidth to avoid hardware limitations.

TR/T1 Effect in Wideband experiments

The proposed Coherent Wideband technique achieves a W-fold acceleration by reducing TR to 1/W of the original TR of a standard scan. Since shortened TR reduces SNR due to the T1 effect, the Wideband technique acquires W times as many points as the standard scan to provide a  $\sqrt{W}$ -fold SNR boost. This process can ensure Wideband results to have comparable SNR as standard results if a sufficiently long TR value is chosen, as the experimental results in this study demonstrate. However, due to the T1 effect, this process won't be as effective for shorter TR values.


To further investigate this phenomenon, the following experiment was conducted to measure and compare the SNR of rat brain images for a standard scan and a Wideband scan when different TR values were chosen. The SNR difference between standard and Wideband scans is depicted as the SNR ratio in the following chart/plot. (As shown in Table.1 & Fig.16)

In 7T environments, T1 is approximately 1900 ms for the gray matter of rat brains. As the results demonstrate, for TR between 3000 and 5000 ms, the Wideband scan and standard scan provide similar SNR values, which is consistent with our rat brain experiments in this study. In contrast, for TR shorter than 3000 ms, the SNR for Wideband scan decreases. In general, for commonly used TR values in rat brain experiments, Wideband MRI provides acceleration benefits and comparable SNR with conventional techniques.

Repetition Time (TR)	Signal to Noise Ratio Between wideband and standard scan (Wideband SNR/Standard SNR)
500 ms	0.7536
1000 ms	0.7996
2000 ms	0.889
3000 ms	0.9726
4000 ms	1.048
5000 ms	1.1151
6000 ms	1.1725
7000 ms	1.2207
8000 ms	1.2606

10000 ms	1.3193
----------	--------

**Table 1.** Signal to noise ratio (SNR) difference between Wideband and standard scans depends on the chosen repetition time (TR).




**Fig. 16. SNR difference between Wideband and standard scans**

Shortened TR reduces SNR due to the T1 effect, the Wideband technique acquires  $W$  times as many points as the standard scan to provide a  $\sqrt{W}$ -fold SNR boost. As the result, the SNR difference between Wideband and standard scans depend on the chosen TR values.

#### The separation limitation of multi-slice distance

The magnitude of separation gradient depends on the distance between two simultaneously-excited slices – proximate slices require stronger gradients for successful separation on the image domain. Since the gradient slew rate is limited, the designed amplitude for additional gradients couldn't be reached for adjacent slices. This results in incorrect phase for acquired signals which shifts the images to undesigned locations and forms artifacts. As shown in Fig. 17, since ME-Wideband (Fig.1) applies 3 times stronger refocusing gradient to offset the accumulated phase, apparent artifacts appeared and became more obvious as slice distance decreased to 1 mm. On the other hand, the Coherent Wideband EPI method utilizes bipolar separation gradients to keep isochromatic spins in phase, which massively reduce the required separation gradient strength. Hence, proximate simultaneously-excited slices can still be acquired and maintained consistent quality.




**Fig. 17. Slice distance effect for ME-Wideband and Coherent Wideband EPI.**

Since ME-Wideband applies stronger refocusing gradients for image separation, its image quality would be limited by the slice distance. As slice distance decreases, ME-Wideband images deteriorate with severe artifacts while Coherent Wideband images maintain consistent quality.

### DTI Precision Analyses

The trade-off between sufficient SNR and image resolution is often a dilemma for MRI researchers. Toward this end, this study presented two usages of Wideband EPI, thin-slice and high-SNR Coherent Wideband DTI, to enhance the calculated MRI image quality and meet the demands for precision medicine studies.

In DTI applications, the partial volume effect causes specious tractography due to fiber kissing or crossing<sup>35</sup>, especially in low-SNR circumstances. As demonstrated in our DTI study results (Fig. 18b), 2X-Thinner Wideband DTI provided robust tracking results, with distinct lateral tracks in the CC region and anteroposterior tracks in the CG region. 2X-Average Wideband DTI also provides smoother tracking results with better reproducibility due to its 1.4x-SNR enhancement (Fig. 18c). On contrary, in the conventional DTI, complex fiber orientations cannot be distinguished with lower spatial resolution, therefore only the dominant lateral tracks remained in the tractography along with some erroneous anteroposterior forks (Fig. 18a). This comparison highlights the importance of Wideband: when the SNR for a DTI scan isn't sufficient to provide a reasonable tractography, Wideband can either enhance SNR or reduce slice thickness to elevate the accuracy of neuro-tractography analyses.



**Fig. 18. Tracking results of conventional and Wideband DTI.**

The tractographies of conventional EPI, 2X-Thinner Wideband and 2X-Average SNR enhanced Wideband

EPI. The SNR of conventional DTI result is 55, slightly lower than previous DTI studies in Fig. 7. In comparison, the 2X-Thinner Wideband result (b) shows better delineation for neuro-architectures, and 2X-Average Wideband result (c) reveals consistent tracking in the corpus callosum (cc) and cingulum (cg) regions.

### Compatibility with current acceleration techniques

The Coherent Wideband EPI technique is highly compatible with multiple existing MRI acceleration techniques, such as Hadamard encoding, and parallel imaging. The Hadamard encoding method adjusts RF phase to reconstruct simultaneously-excited slices, so even more slices can be simultaneously acquired with its combination with Wideband. With multi-coil arrays, Coherent Wideband can also combine with other SMS techniques to reconstruct images with coil sensitivity information and utilize parallel imaging techniques like SENSE and GRAPPA to achieve further acceleration.

## CONCLUSION

This study proposed the Coherent Wideband EPI technique and demonstrated its multifaceted benefits without the usage of phase array coils. This technique also adopted a precise phase alignment method to minimize the  $N/2$  ghost artifacts in conventional EPI images, reducing the ghost factor from 6.3% to 3.1%. The Coherent Wideband EPI technique was employed on a pre-clinical animal system for various applications to demonstrate its versatility. First, a phantom experiment was conducted to verify the consistency between conventional EPI and Wideband EPI images. Subsequently, the structural images for rodent brains with double imaging speed and double spatial resolution were successfully acquired to illustrate how Wideband elevates MRI efficiency. For DTI, under the same scan time, Coherent Wideband EPI provided more accurate and consistent tracking results by either increasing the z-resolution or improving image SNR. Quantitative diffusion analysis (HWHM of FA difference and DA) of the SNR-enhanced Wideband provided the lower fluctuations (FA difference: 56.9%, DA: 64%) for outcome compared to conventional scans.

Regarding the resting fMRI study, the tSNR map was elevated by 44.2% from  $43.28 \pm 9.65$  to  $62.41 \pm 9.52$  by shortening the sampling interval and was boosted by 65.2% to  $70.9 \pm 9.22$  by doubling the scan number in the cortex area. This ends up as a 38.3% increase of associated activation regions in comparison with traditional resting fMRI results due to the tSNR enhancement in the double-sampling case.

In conclusion, the proposed Coherent Wideband EPI technique provides a novel and reliable alternative to achieve simultaneous multi-slice imaging without phase array coils, which could be especially advantageous for small animal pre-clinical MRI scanners. The comprehensive enhancements of this technique not only benefit neuro-architecture imaging, but also advance DTI and fMRI connectivity analysis studies which emphasizes the significance of precise quantitative analyses. Since the thermal noise of MR images may exacerbate error propagation with each step of the DTI/fMRI neuro-estimations, the Wideband technique, which provides credible anatomical information/precision analysis, can greatly contribute to the biological plausibility of neuro-informatics. In summary, Coherent Wideband EPI will provide faster, higher resolution, thinner slice imaging, or higher signal-to-noise MR imaging for precision neuro-architecture studies to facilitate future medical researches.

## References

- 1 Lauterbur, P. C. Image formation by induced local interactions: examples employing nuclear magnetic resonance. *nature* **242**, 190-191 (1973).
- 2 Kwong, K. K. *et al.* Dynamic magnetic resonance imaging of human brain activity during primary sensory stimulation. *Proceedings of the National Academy of Sciences of the United States of America* **89**, 5675-5679, doi:10.1073/pnas.89.12.5675 (1992).
- 3 Bandettini, P. A., Wong, E. C., Hinks, R. S., Tikofsky, R. S. & Hyde, J. S. Time course EPI of human brain function during task activation. *Magnetic Resonance in Medicine* **25**, 390-397, doi:10.1002/mrm.1910250220 (1992).
- 4 Basser, P. J., Mattiello, J. & LeBihan, D. Estimation of the Effective Self-Diffusion Tensor from the NMR Spin Echo. *Journal of Magnetic Resonance, Series B* **103**, 247-254, doi:<https://doi.org/10.1006/jmrb.1994.1037> (1994).
- 5 Basser, P. J., Mattiello, J. & LeBihan, D. MR diffusion tensor spectroscopy and imaging. *Biophys J* **66**, 259-267, doi:10.1016/S0006-3495(94)80775-1 (1994).
- 6 Pierpaoli, C. & Basser, P. J. Toward a quantitative assessment of diffusion anisotropy. *Magnetic Resonance in Medicine* **36**, 893-906, doi:10.1002/mrm.1910360612 (1996).
- 7 Müller, S. Multifrequency selective rf pulses for multislice MR imaging. *Magnetic Resonance in Medicine* **6**, 364-371, doi:10.1002/mrm.1910060315 (1988).
- 8 Souza, S. P., Szumowski, J., Dumoulin, C. L., Plewes, D. P. & Glover, G. SIMA: simultaneous multislice acquisition of MR images by Hadamard-encoded excitation. *J Comput Assist Tomogr* **12**, 1026-1030 (1988).
- 9 Glover, G. H. Phase-offset multiplanar (POMP) volume imaging: a new technique. *J Magn Reson Imaging* **1**, 457-461, doi:10.1002/jmri.1880010410 (1991).
- 10 Larkman, D. J. *et al.* Use of multicoil arrays for separation of signal from multiple slices simultaneously excited. *Journal of Magnetic Resonance Imaging* **13**, 313-317, doi:10.1002/1522-2586(200102)13:2<313::AID-JMRI1045>3.0.CO;2-W (2001).
- 11 Breuer, F. A. *et al.* Controlled aliasing in parallel imaging results in higher acceleration (CAIPIRINHA) for multi-slice imaging. *Magnetic Resonance in Medicine* **53**, 684-691, doi:10.1002/mrm.20401 (2005).
- 12 Setsompop, K. *et al.* Blipped-controlled aliasing in parallel imaging for simultaneous multislice echo planar imaging with reduced g-factor penalty. *Magnetic Resonance in Medicine* **67**, 1210-1224, doi:10.1002/mrm.23097 (2012).
- 13 Xu, J. *et al.* in *Proceedings of the 20th Annual Meeting of ISMRM*.
- 14 Feinberg, D. A. *et al.* Multiplexed echo planar imaging for sub-second whole brain fMRI and fast diffusion imaging. *PLoS One* **5**, e15710-e15710, doi:10.1371/journal.pone.0015710 (2010).
- 15 Moeller, S. *et al.* Multiband multislice GE-EPI at 7 tesla, with 16-fold acceleration using partial parallel imaging with application to high spatial and temporal whole-brain fMRI. *Magnetic Resonance in Medicine* **63**, 1144-1153, doi:10.1002/mrm.22361 (2010).
- 16 Yuen, N. H., Osachoff, N. & Chen, J. J. Intrinsic Frequencies of the Resting-State fMRI Signal: The Frequency Dependence of Functional Connectivity and the Effect of Mode Mixing. *Frontiers in Neuroscience* **13**, doi:10.3389/fnins.2019.00900 (2019).
- 17 Feinberg, D. A., Vu, A. T. & Beckett, A. Pushing the limits of ultra-high resolution human brain imaging with SMS-EPI demonstrated for columnar level fMRI. *NeuroImage* **164**, 155-163, doi:<https://doi.org/10.1016/j.neuroimage.2017.02.020> (2018).
- 18 Hsu, Y.-C. *et al.* Simultaneous multi-slice inverse imaging of the human brain. *Sci Rep* **7**, 17019-17019, doi:10.1038/s41598-017-16976-0 (2017).

- 19 Weaver, J. B. Simultaneous multislice acquisition of MR images. *Magnetic Resonance in Medicine* **8**, 275-284, doi:10.1002/mrm.1910080305 (1988).
- 20 Wu, E. L., Chiueh, T. D. & Chen, J. H. Multiple-frequency excitation wideband MRI (ME-WMRI). *Med Phys* **41**, 092304, doi:10.1118/1.4893502 (2014).
- 21 Wu, E. L. *et al.* Ultra-fast brain MR imaging using simultaneous multi-slice acquisition (SMA) technique. *Annual International Conference of the IEEE Engineering in Medicine and Biology Society. IEEE Engineering in Medicine and Biology Society. Annual International Conference* **2007**, 2618-2621, doi:10.1109/iembs.2007.4352866 (2007).
- 22 Wu, E., Chen, J. & Chiueh, T. in *2010 Annual International Conference of the IEEE Engineering in Medicine and Biology*. 5681-5684 (IEEE).
- 23 Setsompop, K. *et al.* Blipped-controlled aliasing in parallel imaging for simultaneous multislice echo planar imaging with reduced g-factor penalty. *Magn Reson Med* **67**, 1210-1224, doi:10.1002/mrm.23097 (2012).
- 24 Bruder, H., Fischer, H., Reinfelder, H. E. & Schmitt, F. Image reconstruction for echo planar imaging with nonequidistant k-space sampling. *Magn Reson Med* **23**, 311-323, doi:10.1002/mrm.1910230211 (1992).
- 25 Ordidge, R. J. & Cresshull, I. D. The correction of transient B0 field shifts following the application of pulsed gradients by phase correction in the time domain. *Journal of Magnetic Resonance (1969)* **69**, 151-155, doi:[https://doi.org/10.1016/0022-2364\(86\)90228-3](https://doi.org/10.1016/0022-2364(86)90228-3) (1986).
- 26 Beaumont, M., Lamalle, L., Segebarth, C. & Barbier, E. L. Improved k-space trajectory measurement with signal shifting. *Magnetic Resonance in Medicine* **58**, 200-205, doi:<https://doi.org/10.1002/mrm.21254> (2007).
- 27 Feldman, H. M., Yeatman, J. D., Lee, E. S., Barde, L. H. F. & Gaman-Bean, S. Diffusion tensor imaging: a review for pediatric researchers and clinicians. *J Dev Behav Pediatr* **31**, 346-356, doi:10.1097/DBP.0b013e3181dcaa8b (2010).
- 28 Alexander, A. L., Lee, J. E., Lazar, M. & Field, A. S. Diffusion tensor imaging of the brain. *Neurotherapeutics* **4**, 316-329, doi:10.1016/j.nurt.2007.05.011 (2007).
- 29 Lin, C. P., Tseng, W. Y., Cheng, H. C. & Chen, J. H. Validation of diffusion tensor magnetic resonance axonal fiber imaging with registered manganese-enhanced optic tracts. *Neuroimage* **14**, 1035-1047, doi:10.1006/nimg.2001.0882 (2001).
- 30 Fodor, B. & Fingscheidt, T. in *Speech Communication; 10. ITG Symposium*. 1-4.
- 31 Ziegler, S., Jakoby, B. W., Braun, H., Paulus, D. H. & Quick, H. H. NEMA image quality phantom measurements and attenuation correction in integrated PET/MR hybrid imaging. *EJNMMI physics* **2**, 18, doi:10.1186/s40658-015-0122-3 (2015).
- 32 Zhou, W., Bovik, A. C., Sheikh, H. R. & Simoncelli, E. P. Image quality assessment: from error visibility to structural similarity. *IEEE Transactions on Image Processing* **13**, 600-612 (2004).
- 33 Paxinos, G. & Watson, C. *The rat brain in stereotaxic coordinates: hard cover edition*. (Elsevier, 2006).
- 34 Hsu, A.-L. *et al.* Physiological Contribution in Spontaneous Oscillations: An Approximate Quality-Assurance Index for Resting-State fMRI Signals. *PLoS One* **11**, e0148393, doi:10.1371/journal.pone.0148393 (2016).
- 35 Vos, S. B., Jones, D. K., Viergever, M. A. & Leemans, A. Partial volume effect as a hidden covariate in DTI analyses. *Neuroimage* **55**, 1566-1576, doi:10.1016/j.neuroimage.2011.01.048 (2011).

## Acknowledgements

The authors thank the interdisciplinary MRI/MRS laboratory and C.-H. Hsieh of Instrumentation Center, National Taiwan University for MRI experiments. I would like to express my special thanks of gratitude to J.-Wong, M.-Cheng and Prof. C.-Wu for their helpful discussions. This study was supported by Ministry of Science and Technology (MOST 109-2314-B-002-270) and National Taiwan University (MOST 108-2314-B-002 -121 -MY2). The funders have/had no role in study design, data collection and analysis, decision to publish or preparation of the manuscript.



# Figures






Figure 1

Pulse sequences of multiple wideband EPI (a) Blipped Wideband EPI Sequence with constant separation gradients. (b) ME Wideband EPI Sequence with refocusing gradient. (c) Coherent Wideband EPI sequence with bipolar refocusing gradient. (d) Phase plots of multiple EPI sequences. (e) Signal Magnitude of various EPI sequences. (f) Zoom in View of Signal Magnitude




**Figure 2**

EPI Phase Alignment. Multiple  $N/2$  ghost correction method is implemented in (a), (b) & (c). The uncorrected  $N/2$  ghost artifact is shown in (a.3). The artifact is partially removed by the linear correction method in (b.3), and thoroughly removed by the zero-phase compensation method in (c.3).



**Figure 3**


Slice Dependent Phase Alignment. Separation gradients induce phase shift in k-space which leads to  $N/2$  artifacts on the images shown in (a). The artifacts were partially corrected by theoretical phase as shown in (b) & (c). By applying the measured phase, precise offset compensation for Wideband images was ensured as shown in (d) & (e).



**Figure 4**


Phantom results for multiple Wideband EPI methods. (a) Standard EPI results (TR/TE = 500/55 ms). Slices A & B (10 mm from center respectively) were separately excited and acquired. (b) Blipped Wideband EPI (TR/TE = 250/55ms) result. Slices A & B were simultaneously excited and acquired. (c) The ME-Wideband EPI partially mitigated the blur and showed visible inhomogeneity. (d) Coherent Wideband EPI results are consistent with conventional results.





**Figure 5**


Conventional and 2-fold Temporal Wideband in vivo rodent brain study The In Vivo rodent brain images were acquired by a 7T animal MRI with (a) Conventional Spin-Echo EPI (b) Two-fold temporal resolution Coherent Spin-Echo EPI. Coherent Wideband EPI is 2 times faster while the image quality maintains the same as standard images.



**Figure 6**

Conventional and 2-fold spatial resolution Wideband rodent brain results. (a) For standard EPI, the image resolution was  $312 \times 312 \mu\text{m}^2$  and TR=1000ms, so the total scan time was 6 m 40 s. (b) 2X Coherent

Wideband acquired double data points under the same scan time, improving the image resolution to  $156 \times 156 \mu\text{m}^2$ .




**Figure 7**


The 2-fold temporal resolution diffusion study. In the Coherent Wideband 2x temporal resolution diffusion experiment, the acquisition time of the sequence is reduced from 43 min 20 s to 21 min 40 s. The signal to noise and FA value of the 2X fast image kept well at the equal quality of the original one.

## 2X Average Number

- (a) Commercial Spin-Echo EPI      (b) 2X Average Number  
Wideband Spin-Echo EPI



*Average Number = 40*  
*SNR=76.1±8.2*  
*Scan Time = 43min20s*



*Average Number = 80*  
*SNR=105.2±5.8*  
*Scan Time = 43min20s*

- (c) Commercial FA Map



*Whole Brain FA=0.224±0.143*  
*Corpus Callosum FA=0.613±0.052*


- (d) 2X Average Number  
Wideband FA Map



*Whole Brain FA=0.231±0.122*  
*Corpus Callosum FA=0.665±0.047*

Figure 8


Diffusion study of 2X-Average Coherent Wideband EPI images. In the 2X-Average Coherent Wideband MR diffusion experiment, the SNR was boosted 1.4 times from 76 to 105, as expected. The acquisition time of the sequence remains at 43min 20s.



**Figure 9**


Reproducibility analysis of FA and axonal fiber orientation The quantitative analysis of DTI result, including (a) fractional anisotropy (FA) difference and (b) deviation angle. Due to the 1.4x SNR gain, Wideband results have lower HWHM in both histograms, indicating better reproducibility. The half width at half maximum (HWHM) of whole brain FA difference distribution dropped from 0.072 to 0.041(56.9%), and the HWHM of fiber deviation angle histogram decreased from 25.3° to 16.2° (64%).






**Figure 10**

Diffusion Tracking of 2X-Average Coherent Wideband EPI. Images are magnified for the comparison between conventional EPI (SNR=76.1) and SNR enhanced Wideband EPI (SNR=105.2) tractography (2x scan numbers; same scan time). The Wideband result shows higher SNR and uniform tracking in the external capsule (ec), corpus callosum (cc), and cingulum (cg).




**Figure 11**

The diffusion tracking of 2X Thinner Wideband results. Conventional EPI DTI images (SNR=55) and thinner Wideband EPI DTI images (2x thinner slices; same scan time). The conventional low SNR DTI scan yields specious tracking results, particularly in the corpus callosum (cc), and cingulum (cg). The comparison shows that the Wideband result presents more details of neuro-architecture and consistent tracking in fimbria (fi), ec, cc, and cingulum (cg).




**Figure 12**

The 2x thinner wideband diffusion study. DTI FA maps (Fig. 12a & 12b) and DTI Color maps (Fig. 12c & 12d) are magnified to compare the conventional and thin-slice Wideband results. Since the partial volume effect is reduced, the precision of quantitative DTI analysis is the uniformity of tracking are enhanced in the Hippocampal Commissure (HC) and corpus callosum (cc).



**Figure 13**

Temporal signal to noise ratio map. In the resting functional MRI experiment, the tSNR map was boosted by 65.2% to  $70.9 \pm 9.22$  by doubling the scan number, and was elevated from  $43.28 \pm 9.65$  to  $62.41 \pm 9.52$  (44.2%) by shorten the sampling interval in the cortex area. In the corpus callosum and Putamen, by doubling the scan number in the cortex area, tSNR increased by 46.5% and 42.1%; by shortening the sampling interval, tSNR increased from  $34.23 \pm 11.28$  and  $32.12 \pm 13.22$  to  $42.17 \pm 9.31$  (23.2%) and  $37.71 \pm 9.28$  (17.4%).



**Figure 14**

Resting-state functional connectivity. Resting-state functional connectivity with bandpass filtered at 0.01 and 1Hz. Group resting-state functional connectivity maps ( $p < 0.005$ ,  $N = 6$ ) using seed regions in the (a) right motor network, (b) right somatosensory network. The 2X Wideband EPI (Scan Number=240) shows increased (1.38 times as conventional scan) regional connectivity, but 2X Wideband EPI (Scan Number=120) reveals less (0.78 times as conventional scan) connections in the right motor network. In the right somatosensory network, 2X Wideband EPI (Scan Number=240) presents 1.34 times relative regions but 0.82 times regions are observed in 2X Wideband EPI (Scan Number=120).






Figure 15

Multiple applications of the Wideband MRI. This study demonstrated multiple benefits of the Coherent Wideband EPI technique, including temporal/spatial resolution enhancement, SNR improvement, and accuracy study.




Figure 16

SNR difference between Wideband and standard scans Shortened TR reduces SNR due to the T1 effect, the Wideband technique acquires  $W$  times as many points as the standard scan to provide a  $\sqrt{W}$ -fold SNR boost. As the result, the SNR difference between Wideband and standard scans depend on the chosen TR values.





Figure 17

Slice distance effect for ME-Wideband and Coherent Wideband EPI. Since ME-Wideband applies stronger refocusing gradients for image separation, its image quality would be limited by the slice distance. As slice distance decreases, ME-Wideband images deteriorate with severe artifacts while Coherent Wideband images maintain consistent quality.



**Figure 18**

Tracking results of conventional and Wideband DTI. The tractographies of conventional EPI, 2X-Thinner Wideband and 2X-Average SNR enhanced Wideband EPI. The SNR of conventional DTI result is 55, slightly lower than previous DTI studies in Fig. 7. In comparison, the 2X-Thinner Wideband result (b) shows better delineation for neuro-architectures, and 2X-Average Wideband result (c) reveals consistent tracking in the corpus callosum (cc) and cingulum (cg) regions.